

# General Warnings and Cautions

The following general warnings and cautions are applicable to Vescent Photonics instruments.

**WARNING**

This instrument is intended for use by qualified personnel who recognize shock hazards or laser hazards and are familiar with safety precautions required to avoid possible injury. Read the instruction manual thoroughly before using to become familiar with the instrument's operations and capabilities.

**CAUTION**

There are no serviceable parts inside the instrument. Work performed by persons not authorized by Vescent Photonics may void the warranty.

**CAUTION**

Although ESD protection is designed into the instrument, operation in a static-free work area is recommended.

**WARNING**

To avoid electrical shock hazard, connect the instrument to properly earth-grounded, 3-prong receptacles only. Failure to observe this precaution can result in severe injury or death.

**WARNING**

Do not clean outside surfaces of any Vescent Photonics products with solvents such as acetone. Front panels on electronics modules may be cleaned with a mild soap and water solution. Do not clean optics modules.

From:

<https://www.vescent.com/manuals/> - **Product Manuals**

Permanent link:

[https://www.vescent.com/manuals/doku.php?id=warnings\\_cautions&rev=1375286647](https://www.vescent.com/manuals/doku.php?id=warnings_cautions&rev=1375286647)

Last update: **2021/08/26 14:26**

