

The control panel features a vertical sidebar on the left with five icons: a red and blue diamond logo, a home icon, a refresh icon, a lock icon, and a gear icon. The main area is divided into several sections:

- System Status:** Two buttons, "SYSTEM OFF" and "Slow Servo", are located at the top.
- OSC CAVITY:** A horizontal bar with a white line in the center. Below it, "PZT Voltage [V]:" is followed by a slider from -2.0 to 100.0.
- DC Bias [V]:** A section with two columns. The first column has a "PZT" button and a value of "0.00". The second column has a "Temperature" button, a "Target [°C]" value of "25.000", and an "Actual [°C]" value of "-0.600".
- PUMP LASERS:** A horizontal bar with a white line in the center. Below it, "OSC Mod. [mA]:" is followed by a slider from 0.0 to 100.0.
- Current and Temperature:** A section with two columns. The first column has an "OSC" button, a "Curr. [mA]" value of "0", and a "Temp [°C]" value of "-0.200". The second column has an "AMP" button, a "Curr. [mA]" value of "0", and a "Temp [°C]" value of "-0.400".